

JPD-900W Blood Pressure Monitor

(Wrist type)



Features

- ❖ Oscillometric method, automatic inflation & measurement.
- ❖ Fuzzy logic & large LCD display.
- ❖ Heartbeat measurement & Irregular heartbeat detection.
- ❖ Time and date display.
- ❖ 90 memory ;
- ❖ KPa/mmHg Unit exchange.
- ❖ WHO Classification Indicator.
- ❖ PMMA faceplate (acrylic glass)
- ❖ Automatically turns off (within 1 minute) to save power.
- ❖ Measuring Range: pressure: 30mmHg~280mmHg; pulse: 40~190 times/Min.
- ❖ Measuring Accuracy: pressure: +/-3mmHg; pulse: +/-5% of reading.
- ❖ Cuff size: 13.5-19.5cm
- ❖ Unit size: 78x78x35mm
- ❖ LCD size: 44x36mm

فلكه دوم شهران - پلاک ۱۲۷ - مجتمع ماهان - ط. ۲ - واحد ۱۳

☎ ۴۴۳۲۵۹۱۳-۵

📠 ۴۴۳۲۸۷۸۵

✉ info@samanteb.ir

🌐 www.samanteb.ir