

JPD-100S Fetal Doppler



Features

- ❖ Delicate appearance, compact design and portable.
- ❖ Simple one-button operation, controlling power and volume as well.
- ❖ Double-channel audio socket design enables it to provide synchronous audio output signals, accessible for headphones or recording devices.
- ❖ Audio signals can be transmitted to another recording device through attached recording cable .
- ❖ High sensitive ultrasound probe can detect movement sounds of fetal cardiovascular wall at the 16th week of gestation.
- ❖ Unique manufacturing technology is adopted to generate clear fetal heart sound by minimizing ultrasound noise.
- ❖ Ultrasonic frequency: 3.0MHZ
- ❖ Battery: 9V

Packaging information

Unit dimension:

G/W(KG):0.30 N/W(KG):0.20 Size:150×200×70mm

Carton dimension: (30PCS/carton)

G/W(KG):10 N/W(KG):9 Size:480×440×400mm

فکله دوش سهران - پلاک ۱۲۷ - مجتمع مهن - ط. ۲ - واحد ۱۳

☎ ۴۴۲۲۵۹۱۲-۵

📠 ۴۴۲۲۸۷۸۵

✉ info@samanteb.ir

🌐 www.samanteb.ir